



For healthcare professionals

Credenhill is the leading specialist distributor of compression hosiery and equipment. We supply hospitals, surgeries and pharmacies across the UK.

Supplying hosiery direct to patients

Daylong supplies a huge range of fashionable compression hosiery by all leading manufacturers online and by mail order. One stop compression hosiery service from a very experienced team direct to the public.

www.daylong.co.uk

For all your compression hosiery requirements visit daylong.co.uk and step out in style and comfort.

Free Prescription Delivery Service

As an enhanced service, we can dispense your prescription for Open Toe Easy-Slide, Magnide and medical support stockings. So for those of you who live in rural areas or would find it more convenient to have your prescription delivered, send us your prescription and we will process and send the item directly back to your home address.

For more information, please call 0115 932 0144 or visit: www.daylong.co.uk/prescriptions

Call Freephone 0800 195 0650
or visit www.daylong.co.uk
or www.credenhill.co.uk
email: sales@daylong.co.uk
fax: 0115 944 0437

Credenhill Limited

10 Cossall Industrial Estate, Ilkeston,
Derbyshire DE7 5UG England



A Range of Application Aids

for Compression Stockings, Tights & Armsleeves



making your life easier

Available on Prescription

easy-slide

another ARION brand

Putting on Open Toe Compression Stockings and Armsleeves every day has never been easier

To work effectively, compression stockings must apply the right amount of pressure to the leg. This can make them difficult to apply. The Easy-Slide range has been designed to specifically solve these problems.



5 easy steps to fitting open toe stockings

- 1 Slide the Easy-Slide over the foot
- 2 Slide the stocking over the Easy-Slide
- 3 Draw the stocking over the heel
- 4 Grasp the loop on the Easy-Slide
- 5 Pull the Easy-Slide out



Product Range

Easy-Slide for Open Toe Stockings			
Size	PIP Code	Product Code	Shoe Size
Small	266-9935	(EO 1102)	2-4 1/2
Medium	266-9943	(EO 1202)	5-7 1/2
Large	266-9950	(EO 1302)	8-11
Extra Large	266-9968	(EO 1352)	> 11

To select size provide shoe size.

Watch application video at:
www.daylong.co.uk/easy-slide

Available on Prescription from February 2014

magnide

NEW! INNOVATIVE MAGNIDE TO REPLACE THE CARAN

The ultimate donning aid for closed and open toe compression stockings and tights

Magnide provides an effective solution for patients who wear compression stockings and tights with closed and open toe. Hosiery can now be donned without difficulty. Magnide is made from a smooth material that reduces friction significantly helping to minimise the physical strain experienced by clinicians, carers and patients.



Product benefits:

- ✓ AVAILABLE ON FP 10
- ✓ Reduction of physical strain for the clinicians, carer and patient
- ✓ Increase in patient's independence and therapy compliance
- ✓ Stockings retain therapeutic graduated compression longer
- ✓ Easy to use
- ✓ Saves time and money
- ✓ Durable
- ✓ You can take it anywhere

Product Range

Magnide for Open and Closed Toe Stockings			
Size	PIP Code	Product Code	Shoe Size
Medium	380-7534	(EO 1513)	<8
Large	380-7526	(EO 1514)	8-11
Extra Large	387-0375	(EO 1515)	>11

Watch application video at:
www.daylong.co.uk/magnide



easy-slide arm

another ARION brand

An effective solution for patients who wear arm sleeves

Easy-Slide Arm has been specifically designed to alleviate the problems experienced when donning arm sleeves. Easy-Slide Arm is made from very smooth material, which has been treated with a special coating. These qualities create a super-smooth aid that reduces friction significantly.



Product benefits:

- ✓ Reduction of physical strain for clinicians and patients
- ✓ Increase in patient's independence and therapy compliance
- ✓ Easy to use
- ✓ Durable
- ✓ Saves time and money
- ✓ You can take it anywhere
- ✓ Custom-made products available

Product Range

Easy-Slide Arm			
Size	PIP Code	Product Code	Circumference of widest part of the arm
Medium	266-9976	(EO 1402)	Up to 38cm
Large	290-0371	(EO 1425)	38-50cm

Watch application video at:
www.daylong.co.uk/easy-slidearm

